

Tego™ Needlefree Access Device / IV Extension Set

en	- DIRECTIONS
zh	- 用法说明
da	- ANVISNINGER
nl	- AANWIJZINGEN
fr	- MODE D'EMPLOI
de	- GEBRAUCHSANWEISUNG
it	- ISTRUZIONI
jp	- 道順
pt	- INSTRUÇÕES
es	- INSTRUCCIONES
sv	- ANVISNINGAR
pt-br	- ORIENTAÇÕES
pl	- WSKAZÓWKI
cs	- INSTRUKCE
lt	- NURODYMAI
sk	- POKYNY
sl	- NAVODILA
tr	- TALİMATLAR

These products are covered by one or more U.S. and foreign patents. DF-3844 Rev 22 (09/20)

www.icumed.com/symbols-glossary

<u>en – English</u>
Tego™ Needlefree Access Device / IV Extension Set

Use aseptic technique. Single patient use only. Do not resterilize. Device is not made with natural rubber latex. Device is not made with DEHP. Sterile, non-pyrogenic fluid path in unopened undamaged package.

INTENDED USE/INTENDED PURPOSE:
Tego is a needle-free capping device, which closes the end of a catheter, creating a mechanically and microbiologically closed system when attached to the hub of a catheter. Tego will permit access to the catheter without the use of needles, and therefore, passively aid in the reduction of needlestick injuries.

INDICATIONS FOR USE:
Tego is intended for use as an accessory to a vascular access device (catheter) used in hemodialysis and apheresis applications or as an accessory to an intravascular administration set for the administration or withdrawal of fluids to a patient through a cannula or needle placed in the vein or artery.

CONTRAINDICATIONS:
None known.

WARNINGS - RESIDUAL RISK:
Inspect device and read Instructions for Use prior to usage. A damaged device or improper use may lead to the following hazardous situations: delay of therapy, blood loss, ergonomic stress, embolism from air/blood clot/particulate, exposure to allergenic substances/biological contaminants/toxins/blood products/patient blood, infusion of drug precipitate. Harms potentially associated with these hazardous situations include reversible injury, permanent injury, or death.

- PRECAUTIONS:**
- Do not use needles, blunt cannulas, or luer caps on Tego.
 - Tego is compatible with ISO 80369-7 standard male luers (Note: specific to item number D1005 and 011-D1005 only)
 - Disinfect Tego before each access using 70% isopropyl alcohol in aggressive circular motion for three (3) seconds or per facility protocol.
 - Access Tegoby pushing straight on into the connector (without angled or sliding entry). Do not twist or bend the connector when accessing.
 - Tighten and secure all non-bonded connections and inspect for leaks before use. If connections are bonded, do not attempt to tighten or disconnect. Doing so may damage or break the connection.
 - Clamp line before disconnecting from Tego, and between dialysis sessions per established clinical and catheter manufacturer’s practices
 - Device should be changed in accordance with current, recognized guidelines of IV therapy.
 - Tego has been validated for functional and microbiological efficacy for up to seven days.

- INSTRUCTIONS FOR USE:**
- Using aseptic technique, remove Tego from package.
 - Prime Tego in accordance with facility procedure/protocol.
 - Remove cap and attach Tego to desired vascular access device.
 - Swab silicone seal of Tego before each use in accordance with facility protocol.
 - Attach administration device or syringe by pushing straight into Tego access device for infusion. **When using a luer slip connection, only insert half-way while turning slightly clockwise. When removing, apply the same clockwise turning motion. Do not overtighten.**
 - Clamp vascular access device before disconnecting a device (blood line, blood withdrawal device, syringe).
 - Dispose of device in accordance with facility protocol.

- CAUTIONS:**
- Do not reuse or resterilize. Reuse negatively impacts performance / sterility potentially resulting in product failure / contamination.
 - Federal (USA) law restricts the use of this device to sale by or on the order of a physician (for USA only).
 - Device is intended to be used and handled by trained clinical personnel only. Care should be exercised to prevent unintended handling by children.
 - Do not use if package is damaged or opened prior to use.
 - Use by expiration date on package label.

NOTES:

- Straight fluid path permits flow rates of greater than 600 mL/min.
- Tego is comprised of non-metal components and is fully compatible with MRI procedures.
- Use of Tego may reduce the risk of catheter-related bloodstream infections (CRBSI) [1].
- Using saline as a lock solution with Tego may reduce the use of heparin or alternative antithrombotic locking solutions [1, 2, 3].

ADVERSE EVENTS:
Use of the Tego is consistent with the current knowledge and state of art for devices used in infusion therapy. In normal use conditions, the clinical benefits of Tego outweigh the associated risks.

Serious incidents that occur during use of the product should be reported to ICU Medical. These reports can be submitted at www.icumed.com. Incidents should also be reported to the relevant local regulatory authority where required.

- BIBLIOGRAPHY:**
- Brunelli SM, Njord L, Hunt AE, Sibbel SP. Use of the Tego needlefree connector is associated with reduced incidence of catheter-related bloodstream infections in hemodialysis patients. Int J Nephrol Renovasc Dis. 2014;7:131-139. Published 2014 Apr 3. doi:10.2147/IJNRD.S59937.
 - Farthing C, Krishnan M, Mayne T, Collard S, Nissenson A. Tego connectors reduce heparin use without affecting blood flow rate compared to traditional central venous catheter locks. American Nephrology Nurses’ Association 43rd National Symposium. Orlando, FL, Apr 29-May 2, 2012.

- Krishnan M, Mayne T, Farthing C, Collard S,Nissenson A. TEGO Connectors Reduce Heparin Use without Affecting Blood Flow Rate Compared to Traditional Central Venous Catheter Locks. American Society of Nephrology Kidney Week. Philadelphia, PA, Nov 8-11, 2011.

These directions for use were drawn up in English and have been translated into other languagues. In the event of inconsistencies of interpretations, the English text shall be determinative.

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BIBLIOGRAFIE:

- Brunelli SM, Njord L, Hunt AE, Sibbel SP. Use of the Tego needlefree connector is associated with reduced incidence of catheter-related bloodstream infections in hemodialysis patients. Int J Nephrol Renovasc Dis. 2014;7:131-139. Gepubliceerd op 3 april 2014. doi:10.2147/IJNRD.S59937.
- Farthing C, Krishnan M, Mayne T, Collard S, Nissenson A. Tego connectors reduce heparin use without affecting blood flow rate compared to traditional central venous catheter locks. American Nephrology Nurses' Association 43rd National Symposium. Orlando, FL, 29 april-2 mei 2012.
- Krishnan M, Mayne T, Farthing C, Collard S,Nissenson A. TEGO Connectors Reduce Heparin Use without Affecting Blood Flow Rate Compared to Traditional Central Venous Catheter Locks. American Society of Nephrology Kidney Week. Philadelphia, PA, 8-11 november 2011.

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BIBLIOGRAPHIE:

- Brunelli SM, Njord L, Hunt AE, Sibbel SP. Use of the Tego needlefree connector is associated with reduced incidence of catheter-related bloodstream infections in hemodialysis patients. Int J Nephrol Renovasc Dis. 2014;7:131-139. Published 2014 Apr 3. doi:10.2147/IJNRD.S59937.
- Farthing C, Krishnan M, Mayne T, Collard S, Nissenson A. Tego connectors reduce heparin use without affecting blood flow rate compared to traditional central venous catheter locks. American Nephrology Nurses' Association 43rd National Symposium. Orlando, FL, Apr 29-May 2, 2012.
- Krishnan M, Mayne T, Farthing C, Collard S,Nissenson A. TEGO Connectors Reduce Heparin Use without Affecting Blood Flow Rate Compared to Traditional Central Venous Catheter Locks. American Society of Nephrology Kidney Week. Philadelphia, PA, Nov 8-11, 2011.

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BIBLIOGRAFIA:

- Brunelli SM, Njord L, Hunt AE, Sibbel SP. Use of the Tego needlefree connector is associated with reduced incidence of catheter-related bloodstream infections in hemodialysis patients. Int J Nephrol Renovasc Dis. 2014;7:131-139. Pubblicato il 3 aprile 2014. doi:10.2147/IJNRD.S59937.
- Farthing C, Krishnan M, Mayne T, Collard S, Nissenson A. Tego connectors reduce heparin use without affecting blood flow rate compared to traditional central venous catheter locks. American Nephrology Nurses' Association 43rd National Symposium. Orlando, FL, 29 aprile-2 maggio 2012.
- Krishnan M, Mayne T, Farthing C, Collard S, Nissenson A. TEGO Connectors Reduce Heparin Use without Affecting Blood Flow Rate Compared to Traditional Central Venous Catheter Locks. American Society of Nephrology Kidney Week. Filadelfia, PA, 8-11 novembre 2011.

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BIBLIOGRAFIA:

- Brunelli SM, Njord L, Hunt AE, Sibbel SP. Use of the Tego needlefree connector is associated with reduced incidence of catheter-related bloodstream infections in hemodialysis patients. Int J Nephrol Renovasc Dis. 2014; 7: 131-139. Publicado em 3 de abril de 2014. doi: 10.2147/IJNRD.S59937.
- Farthing C, Krishnan M, Mayne T, Collard S, Nissenson A. Tego connectors reduce heparin use without affecting blood flow rate compared to traditional central venous catheter locks. 43.º Simpósio Nacional da Associação Americana de Enfermeiros em Nefrologia. Orlando, FL, 29 de abril a 2 de maio de 2012.
- Krishnan M, Mayne T, Farthing C, Collard S, Nissenson A. TEGO Connectors Reduce Heparin Use without Affecting Blood Flow Rate Compared to Traditional Central Venous Catheter Locks. Semana do Rim da Sociedade Americana de Nefrologia. Filadélfia, PA, 8 a 11 de novembro de 2011.

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LITERATURA:

- Brunelli SM, Njord L, Hunt AE, Sibbel SP. Use of the Tego needlefree connector is associated with reduced incidence of catheter-related bloodstream infections in hemodialysis patients. Int J Nephrol Renovasc Dis. 2014;7:131-139. Objavljeno 3. aprila 2014. doi:10.2147/IJNRD.S59937.
- Farthing C, Krishnan M, Mayne T, Collard S, Nissenson A. Tego connectors reduce heparin use without affecting blood flow rate compared to traditional central venous catheter locks. American Nephrology Nurses’ Association 43rd National Symposium. Orlando, FL, 29. april–2. maj 2012.
- Krishnan M, Mayne T, Farthing C, Collard S,Nissenson A. TEGO Connectors Reduce Heparin Use without Affecting Blood Flow Rate Compared to Traditional Central Venous Catheter Locks. American Society of Nephrology Kidney Week. Philadelphia, PA, 8. –11. november 2011.

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KAYNAKÇA:

- Brunelli SM, Njord L, Hunt AE, Sibbel SP. Use of the Tego needlefree connector is associated with reduced incidence of catheter-related bloodstream infections in hemodialysis patients. Int J Nephrol Renovasc Dis. 2014; 7: 131-139. Published 2014 Apr 3. doi:10.2147/IJNRD.S59937.
- Farthing C, Krishnan M, Mayne T, Collard S, Nissenson A. Tego connectors reduce heparin use without affecting blood flow rate compared to traditional central venous catheter locks. American Nephrology Nurses’ Association 43rd National Symposium. Orlando, FL, Apr 29-May 2, 2012.
- Krishnan M, Mayne T, Farthing C, Collard S,Nissenson A. TEGO Connectors Reduce Heparin Use without Affecting Blood Flow Rate Compared to Traditional Central Venous Catheter Locks. American Society of Nephrology Kidney Week. Philadelphia, PA, Nov 8-11, 2011.

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