



MicroClave™ Clear

Neutral Displacement Connector



Disinfect to protect

- › When placing a new MicroClave Clear connector onto a catheter, open the package using aseptic technique.
- › Attach the MicroClave Clear connector onto the catheter by twisting until tightened.
- › Disinfect the MicroClave Clear connector in accordance with facility protocol for appropriate scrubbing and dry times.



Administer or aspirate

- › Attach the IV tubing, syringe, or blood tube holder to the MicroClave Clear connector by inserting the luer and twisting until a friction fit is achieved.
- › Do not overtighten a luer beyond the friction fit as this may damage both the luer and the MicroClave Clear connector.



Flush after each use

- › Flush the MicroClave Clear connector with normal saline or in accordance with facility protocol. After blood use, the MicroClave Clear connector can be flushed clean and does not require change-out.
- › Use routine flushing in accordance with facility protocol in order to maintain catheter patency.
- › Change the MicroClave Clear connector in accordance with facility protocol and CDC guidelines.

Functional attributes

- › Lipid and blood compatible
- › Radiographic imaging compatible (CT compatible)
- › Contains no latex, phthalates (DEHP) or metal components MRI Compatible

Tips

- › To disconnect, grasp the MicroClave Clear connector and then twist the mating luer away from the MicroClave Clear connector until loose.
- › Do not hold the catheter hub during disconnect as this may cause accidental removal of the MicroClave Clear connector from the hub.



Also available in blue



ANTT® Evaluated Technology promotes safe and effective Aseptic Non-Touch Technique principles by design

Warning: Clave™ connectors may be incompatible with some male-luer connectors including prefilled glass syringes. To avoid damage to the Clave or syringes or male luers which may result in delays of medication administration and possible serious adverse events, users should confirm mating luers or syringes have an internal diameter range of 0.062" to 0.110". Check the internal diameter of the male-luer connector of the mating syringe prior to using it to access the Clave. Products outside of these dimensional tolerances should not be used.



MicroClave Clear connector's saline flush option is designed to help reduce the risk of heparin-induced thrombocytopenia.

icumedical
human connections