

ICU Medical MedNet™ Safety Software
IHE Integration Statement

Vendor	Product Name	Version	Date
ICU Medical	ICU Medical Mednet™	6.32	January 2023

This product implements all transactions required in the IHE technical framework to support the IHE integration profiles, actors, and options listed below:

Integrated Profiles Implemented	Actors Implemented
CT – Consistent Time	Time Client
DEC – Device Enterprise Communication	Device Observation Reporter
PIV – Point-of-Care Infusion Verification	Infusion Order Consumer
ACM – Alert Communication Management	Alert Reporter
IPEC – Infusion Pump Event Communication	Device Observation Reporter
MEMDMC – Medical Equipment Management Device Management Communication	Device Management Information Reporter
Point-of-Care Infusion Verification	Infusion Order Consumer

Internet Address for Vendor’s IHE Information
https://www.icumed.com/products/infusion-therapy/infusion-pumps-and-software/infusion-safety-software/icu-medical-mednet-safety-software/

Links to Standards Conformance Statements for Implementation	
HL7	Not Applicable
DICOM	Not Applicable

Links to General Information on IHE	
In North America:	In Europe:
www.ihe.net	www.ihe-europe.net

