

Warmer Patients.

Level 1™ Snuggle Warm™

Patient Warming Blankets Reference Cards



icumedical
human connections


Snuggle Warm Convective Warming Blankets Reference Cards

At ICU Medical, we believe that innovation goes well beyond products. That is why we collaborate closely with clinicians to ensure that all of our Level 1 products and services satisfy our customers' needs. Level 1 convective warming products provide a comprehensive selection of innovative blanket styles for a wide variety of surgical procedures. Throughout the perioperative period, Snuggle Warm blankets are proven and dependable.

Snuggle Warm Convective Warming Blankets are:

- › Made of soft fabric-like material that is comfortable next to patients' skin — no hot plastic
- › Strong and durable — resist tears, punctures and fluids
- › Designed to conform to the patient, reducing the need for tie-downs
- › Made of foldable blanket material to allow for easy access to the patient
- › Equipped with precisely positioned air holes and channel design to maximize heat delivery
- › Latex-free and radiotranslucent
- › Available in multiple configurations for a broad set of clinical applications

ADULTS BLANKETS




SWU-2001
Overbody
Full Body




SWU-2003
Overbody
Upper Body




SWU-2004
Overbody
Lower Body




SWU-2007
Tube Blanket



SWU-2010
Overbody Small
Upper Body



SWU-2013
Underbody
Full Body



SWU-2016
Overbody
Full Body
Split Access



SWU-2113
Underbody
Full Body with
Arm Openings



SWU-2119
Underbody Poncho
with Arm Openings

PEDIATRIC BLANKETS



SWU-2002
Over Body
Full Body




SWU-2009
Underbody
Full Body



SWU-2011
Underbody
Full Body

Content

SWU-2001	Overbody – Full Body Blanket	1
SWU-2003	Overbody – Upper Body Blanket	2
SWU-2004	Overbody – Lower Body Blanket	3
SWU-2007	Tube Blanket	4
SWU-2010	Over body – Small Upper Body Blanket	5
SWU-2013	Underbody – Full Body Blanket	6
SWU-2016	Overbody – Full Body Split Access Blanket	7
SWU-2113	Underbody – Full Body with Arm Opening Blanket	8
SWU-2119	Underbody – Poncho with Arm Openings Blanket	9
SWU-2002	Overbody – Full Body Pediatric Blanket	10
SWU-2009	Underbody – Full Body Pediatric Blanket	11
SWU-2011	Underbody – Full Body Pediatric Blanket	12




SWU-2001 | Overbody - Upper Body Blanket

- › 102 x 203 cm
- › 80 x 40 in
- › Includes non-inflating foot drape

Key attribute:

- › Large size ensures maximum coverage for all patient sizes




SWU-2003 | Overbody - Upper Body Blanket

- › 203 x 102 cm
- › 80 x 40 in
- › Includes non-inflating flaps

Key attributes:

- › Accordion-style chest panel accommodates large and small patients
- › Clear plastic drape captures and delivers warm air around the patient's head




SWU-2004 | Overbody – Lower body blanket

- › 165 x 102 cm
- › 65 x 40 in
- › Includes non-inflating foot drape

Key attribute:

- › Non-inflating foot drape protects patient's feet




SWU-2007 | Tube Blanket

- › 77 x 178 cm
- › 30 x 70 in

Key attribute:

- › Allows for full access to the patient




SWU-2010 | Overbody – Small Upper Body Blanket

- › 198 x 78 cm
- › 78 x 31 in

Key attributes:

- › Tie-downs aid in maintaining the blanket's position
- › Smaller size allows for increased access to the patient
- › Clear plastic drape captures and delivers warm air around the patient's head




SWU-2013 | Underbody - Full Body Blanket

- › 102 x 203 cm
- › 40 x 80 in

Key attributes:

- › “Halo” design headrest area allows for warming of patient’s head
- › Underbody blanket allows for full access to the patient
- › Three hose ports for convective warmer placement flexibility




SWU-2016 | Overbody - Full Body Split Access Blanket

- › 102 x 203 cm
- › 40 x 80 in

Key attribute:

- › Perforation up the middle of the blanket allows for increased patient access




SWU-2113 | Underbody - Full Body with Arm Openings Blanket

- > 102 x 203 cm
- > 40 x 80 in

Key attributes:

- > Pre-opened arm and head openings save time
- > Wide arm openings allow for easy positioning with arm(s) out or toboggan-style arm restraints
- > Easily accommodates patients in a supine, lateral or prone position
- > Tuck flaps secure the blanket to prevent sliding when transferring the patient to the operating table
- > Underbody blanket allows for full access to the patient
- > User-positioned head drape captures and delivers warm air around the patient's head and neck




SWU-2119 | Underbody - Poncho with Arm Openings Blanket

- › 102 x 164 cm
- › 40 x 65 in

Key attributes:

- › Reversible design allows use pre-operatively and postoperatively as an overbody blanket and during surgery as an underbody blanket
- › Underbody use during surgery enables simultaneous warming of the patient's backside and the front of the upper torso
- › Wide arm openings permit easy positioning with arm(s) out or toboggan-style arm restraints




SWU-2002 | Overbody - Full Body Pediatric Blanket

- › 102 x 146 cm
- › 40 x 58 in
- › Includes non-inflating flaps

Key attribute:

- › Sized to fit pediatric patients




SWU-2009 | Underbody - Full Body Pediatric Blanket

- › 67 x 94 cm
- › 26 x 37 in
- › Includes non-inflating flaps

Key attributes:

- › Underbody blanket allows for full access to the patient
- › Two user-positioned clear plastic drapes to capture warmth in specific locations while allowing access where necessary
- › For patients up to 9 kilograms
- › Uninflated headrest for better access
- › Designed with absorbent pads to minimise pooling of fluid on the blanket
- › Can be used on an operating table, bed or bassinot




SWU-2011 | Underbody - Full Body Pediatric Blanket

- › 102 x 137 cm
- › 40 x 54 in
- › Includes non-inflating flaps

Key attributes:

- › Underbody blanket allows for full access to the patient
- › Two user-positioned clear plastic drapes to capture warmth in specific locations while allowing access where necessary
- › For patients up to 36 kilograms
- › Designed with absorbent pads to minimise pooling of fluid on the blanket
- › Can be used on an operating table, bed or bassinet

ADULTS BLANKETS




SWU-2001
Overbody
Full Body




SWU-2003
Overbody
Upper Body




SWU-2004
Overbody
Lower Body




SWU-2007
Tube Blanket



SWU-2010
Overbody Small
Upper Body



SWU-2013
Underbody
Full Body



SWU-2016
Overbody
Full Body
Split Access



SWU-2113
Underbody
Full Body with
Arm Openings



SWU-2119
Underbody Poncho
with Arm Openings

PEDIATRIC BLANKETS



SWU-2002
Over Body
Full Body



SWU-2009
Underbody
Full Body



SWU-2011
Underbody
Full Body

Rx Only For safe and proper use, refer to the Instructions for Use.

© 2024 ICU Medical, Inc. | P24-6872

icumedical
human connections

800.824.7890 | www.icumed.com