

# PharmGuard™ Infusion Management System

## IHE Integration Statement

Vendor	Product Name	Version	Date
ICU Medical	PharmGuard	2.7	June 2024

This product implements all transactions required in the IHE technical framework to support the IHE integration profiles, actors, and options listed below:

Integrated Profiles Implemented	Actors Implemented
ACM – Alert Communication Management	Alert Reporter
CT – Consistent Time	Time Client
DEC – Device Enterprise Communication	Device Observation Reporter
IPEC – Infusion Pump Event Communication	Device Observation Reporter
MEMDMC – Medical Equipment Management Device Management Communication	Device Management Information Reporter
MEMLS – Medical Equipment Management Location Services	Location Observation Consumer
PIV – Point-of-Care Infusion Verification	Infusion Order Consumer
Point-of-Care Infusion Verification	Infusion Order Consumer

### Internet Address for Vendor’s IHE Information

<https://www.icumed.com/products/infusion-therapy/infusion-pumps-and-software/infusion-safety-software/pharmguard-infusion-management-system/>

### Links to Standards Conformance Statements for Implementation

HL7	Not Applicable
DICOM	Not Applicable

### Links to General Information on IHE

In North America:	In Europe:
<a href="http://www.ihe.net">www.ihe.net</a>	<a href="http://www.ihe-europe.net">www.ihe-europe.net</a>