



Important

Please open this form in Adobe Acrobat.
It has interactive elements which may not function correctly
if opened in your web browser.

If you don't have Adobe Acrobat Reader installed, you can download it here:

<https://get.adobe.com/reader/>

Bivona™ Customized FlexTend™ Tracheostomy Tube Template

1. Shaft Length



- 1a. Shaft Length _____ mm (5 mm min, 60 mm max)
 1b. Shaft Length _____ mm (30 mm min, 140 mm max)

2. Shaft

I.D.	O.D.	I.D.	O.D.	I.D.	O.D.	I.D.	O.D.	I.D.	O.D.
<input type="radio"/> 2.5	4.0	<input type="radio"/> 4.0	6.0	<input type="radio"/> 5.5	8.0	<input type="radio"/> 7.0	10.0	<input type="radio"/> 8.5	11.8
<input type="radio"/> 3.0	4.7	<input type="radio"/> 4.5	6.7	<input type="radio"/> 6.0	8.7	<input type="radio"/> 7.5	10.4	<input type="radio"/> 9.0	12.3
<input type="radio"/> 3.5	5.3	<input type="radio"/> 5.0	7.3	<input type="radio"/> 6.5	9.4	<input type="radio"/> 8.0	11.0	<input type="radio"/> 9.5	13.3


3. Connector (check one box only)

- 3a. Swivel 2.5–9.5 mm 3b. Fixed 2.5–6.0 mm



4. Neck Flange (check one box only)


- 4a. 2.5–6.0 mm 4c. 2.5–9.5 mm 4d. 6.5–9.5 mm 4e. 2.5–4.0 mm 4f. 2.5–6.0 mm



Neck Flange Length: 60 mm Neck Flange Length: 83 mm Neck Flange Length: 67 mm Neck Flange Length: 47 mm Neck Flange Length: 60 mm

5. Cuff (check one box only)


- 5a. TTS Tight to Shaft 5b. Aire-Cuf™ Pedi 5b. Aire-Cuf Adult 5c. Fome-Cuf™ Adult 5d. Cuffless




2.5–9.5 mm 2.5–5.5 mm 5.0–9.5 mm 2.5–9.5 mm 2.5–9.5 mm

6. Shaft Option

- 6a. Shaft (Only proximal section is wire reinforced.)



- 6b. Hyperflex



2.5–9.5 mm

7. Optional

- 7a. Proximal (Top) Cuff

- TTS Cuff
 Aire-Cuf
 Fome-Cuf
 Pediatric
 Adult

Cuff position from Neck Flange _____ mm
 Note: Minimum distance 2 mm

- 7b. Stoma Seal

- Solid Silicone
 _____ mm O.D.
 _____ mm Length
 Tapered
 Non Tapered

- 7c. Distal (Bottom) Cuff

Cuff distance from tip to shaft _____ mm
 Note: Cuff distance from tip cannot be less than the minimum standard distance listed below:
 a. Pediatric TTS, Aire-Cuf & Fome-Cuf – 3 mm
 b. Adult Aire-Cuf & Fome-Cuf – 5 mm
 c. Adult TTS – 9 mm

- 7d. Talk Attachment
 5.0 mm I.D. Adult shaft or above only

8. Approval Please note: ICU Medical will review the completed template for compliance.

Any required modifications will be made with the approval of a clinician.

Hospital/Clinic _____
 Clinician Name _____ Clinician Title _____
 P.O. Number _____ Date _____

9. Select Type of Service

- Extended next day service (non-sterile), same day shipping if order received by 12 noon (CST) USA only
 Standard service (sterile) Express service (non-sterile)
 Standard service (non-sterile) Next day service (non-sterile) USA only

icumedical
 human connections
 Phone: 1-800-424-8662
 Fax: 219-989-7435
 Email: bco@icumed.com
 Web: www.icumed.com