


Gripper Plus™ Safety Port Access Needle



Indications

This product is indicated for the administration into, or withdrawal of fluids from implanted ports. It is designed to help protect against exposure to bloodborne pathogens caused by accidental needlestick injuries. It does not protect against other routes of bloodborne pathogen transmission.

Legal Manufacturer

ICU Medical, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442 USA

Manufacturing Site Name and Address

ICU Medical Healthcare Manufacturing S.A. de C.V.
Avenida Calidad No. 4
Parque Industrial Internacional Tijuana
Tijuana, B.C.C.P. Baja CA 22425 Mexico

Country of Origin

Mexico

Animal Origin

None of the devices listed on this product family incorporate any substances that are derived from animal tissues as per Commission Regulation (EU) 722/2012.

Sterilisation Method

Ethylene Oxide


Shelf-Life

A 5-year expiration date is assigned to the product as long as the packaging is undamaged and unopened, based on accelerated aging stability shelf-life testing and real-time aging.


Not manufactured with DEHP plasticiser

MR conditional up to 3 Tesla


Priming Volume



Without Y-Site Straight Tubing



Y-Site Tubing



Needleless Y-Site Tubing

Labelling and Packaging

| Container Type | Length | Width | Height |
|--------------------|--------|--------|--------|
| Unit Pack (1 unit) | 101 mm | 145 mm | 25 mm |
| Shelf Pack | 205 mm | 92 mm | 148 mm |
| Case Pack | 492 mm | 422 mm | 406 mm |

Component Compositions

| | |
|---------------------|-------------------|
| Clip | Acetal |
| Foam | EVA |
| Base & Pivot Arm | Acrylic Copolymer |
| Cannula | 304 SS |
| Cannula Protector | LDPE |
| Tubing | PVC, TOTM |
| Clamp | Acetal |
| Female Luer | PVC, TOTM |
| Non-Vented Male Cap | ABS |

Product Range and Specifications

Gripper Plus Safety Port Access Needles

Not made with natural rubber latex, non-coring, non-siliconised, removable contoured grip, cushioned needle platform, 20cm (8 in.) extension tubing with or without Y-site

| Needleless Y-site – Luer-activated needleless injection site | Length | Reorder No. | GTIN No. | Unit |
|---|------------------|-------------|-------------------|--------|
| 22 gauge | 16 mm (5/8 inch) | 21-2960-24 | 3 0610586 02768 4 | 12/box |
| | 19 mm (3/4 inch) | 21-2961-24 | 3 0610586 02769 1 | 12/box |
| | 25 mm (1 inch) | 21-2962-24 | 3 0610586 02770 7 | 12/box |
| | 32 mm (1 ¼ inch) | 21-2963-24 | 3 0610586 02771 4 | 12/box |
| 20 gauge | 16 mm (5/8 inch) | 21-2964-24 | 3 0610586 02772 1 | 12/box |
| | 19 mm (3/4 inch) | 21-2965-24 | 3 0610586 02773 8 | 12/box |
| | 25 mm (1 inch) | 21-2966-24 | 3 0610586 02774 5 | 12/box |
| | 32 mm (1 ¼ inch) | 21-2967-24 | 3 0610586 02775 2 | 12/box |
| 19 gauge | 19 mm (3/4 inch) | 21-2968-24 | 3 0610586 02776 9 | 12/box |
| | 25 mm (1 inch) | 21-2969-24 | 3 0610586 02777 6 | 12/box |
| | 32 mm (1 ¼ inch) | 21-2970-24 | 3 0610586 02778 3 | 12/box |
| Y-site – Removable split septum injection cap on Luer-lock Y-site | Length | Reorder No. | GTIN No. | Unit |
| 22 gauge | 16 mm (5/8 inch) | 21-2860-24 | 3 0610586 02564 2 | 12/box |
| | 19 mm (3/4 inch) | 21-2861-24 | 3 0610586 02565 9 | 12/box |
| | 25 mm (1 inch) | 21-2862-24 | 3 0610586 02566 6 | 12/box |
| | 32 mm (1 ¼ inch) | 21-2863-24 | 3 0610586 02567 3 | 12/box |
| 20 gauge | 16 mm (5/8 inch) | 21-2864-24 | 3 0610586 02568 0 | 12/box |
| | 19 mm (3/4 inch) | 21-2865-24 | 3 0610586 02569 7 | 12/box |
| | 25 mm (1 inch) | 21-2866-24 | 3 0610586 02570 3 | 12/box |
| | 32 mm (1 ¼ inch) | 21-2867-24 | 3 0610586 02571 0 | 12/box |
| 19 gauge | 19 mm (3/4 inch) | 21-2868-24 | 3 0610586 02572 7 | 12/box |
| | 25 mm (1 inch) | 21-2869-24 | 3 0610586 02573 4 | 12/box |
| | 32 mm (1 ¼ inch) | 21-2870-24 | 3 0610586 02574 1 | 12/box |
| Without Y-site | Length | Reorder No. | GTIN No. | Unit |
| 22 gauge | 16 mm (5/8 inch) | 21-2770-24 | 3 0610586 02563 5 | 12/box |
| | 19 mm (3/4 inch) | 21-2766-24 | 3 0610586 02559 8 | 12/box |
| | 25 mm (1 inch) | 21-2760-24 | 3 0610586 02553 6 | 12/box |
| | 32 mm (1 ¼ inch) | 21-2761-24 | 3 0610586 02554 3 | 12/box |
| 20 gauge | 16 mm (5/8 inch) | 21-2769-24 | 3 0610586 02562 8 | 12/box |
| | 19 mm (3/4 inch) | 21-2767-24 | 3 0610586 02560 4 | 12/box |
| | 25 mm (1 inch) | 21-2762-24 | 3 0610586 02555 0 | 12/box |
| | 32 mm (1 ¼ inch) | 21-2763-24 | 3 0610586 02556 7 | 12/box |
| 19 gauge | 19 mm (3/4 inch) | 21-2768-24 | 3 0610586 02561 1 | 12/box |
| | 25 mm (1 inch) | 21-2764-24 | 3 0610586 02557 4 | 12/box |
| | 32 mm (1 ¼ inch) | 21-2765-24 | 3 0610586 02558 1 | 12/box |

For more information, visit www.icumed.com or call 866.488.6088

